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Temperatur des Meßgutes: .....

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Elektroden-Einsatzdauer: .....

## Instructions for the use of measurement and reference electrodes

### Information on use, treatment and warranty!

## B 20.2900

### Note

If any difficulties should arise during placing the electrodes into operation, please do not carry out any manipulations on the electrode which are not permitted. You could endanger your rights under the warranty. Please contact the nearest office or the main factory. Any JUMO electrodes which has come to the end of its life and is returned will be disposed of by us in accordance with the appropriate regulations.

The following details are required for dealing with complaints: (Please submit copy!)

Serial number: .....

Electrode cannot be calibrated

Electrode drifts

Other faults: .....

Solution in which electrode is used: .....

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### I. General Notes

1. Electrode systems for determining pH or redox potential consist of a measurement and a reference electrode. Combination pH and redox electrodes contain a measurement and a reference electrode in a common shaft.
2. All electrodes are supplied fitted with a soaking cap. The cap contains de-ionised water (glass electrodes) or 3 M KCl solution (reference electrodes and combination electrodes) in order to ensure that the electrode is always ready for use. The soaking cap has to be removed before use.
3. With regard to plug connections it is essential to ensure that they are kept absolutely clean and dry in order to avoid electrical creepage currents. When the user carries out any installation work on low-noise coaxial cables it is important to note that the black layer between the screening braiding and the inner insulation is semi-conducting and has to be removed completely.
4. When using electrodes with plastic shaft the user has to check that the proposed cleaning method is compatible with the material.

### II. Notes on Electrodes

#### 1. Glass electrodes

- 1.1 The pH-sensitive membrane glass has to be treated with care and protected against damage.
- 1.2 An important condition for satisfactory operation of the glass electrode is the presence of a water-containing layer, a so-called soak layer, on the surface of the glass membrane. If the electrode has been stored dry for a longer period it is necessary to condition it before use. For this purpose it is immersed for about 24 hours in a 3 M KCl solution.
- 1.3 The internal buffer has to cover the internal surface of the glass membrane. Air bubbles in the membrane chamber have to be removed by gently shaking down the electrode in the vertical position (similar to a clinical thermometer). Electrodes have to be installed vertically downwards. The angle to the vertical must not exceed 80°.

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- 1.4 Any contamination deposited on the surface of the glass membrane has to be removed. If careful dabbing with a moist soft paper tissue is not successful, various chemical cleaning methods may be used depending on the type of contamination (gentle glass cleaning agents, laboratory detergents, acetone, alcohol, acid solutions of low strength such as 10% hydrochloric acid). Special cleaning solutions are available for removing protein-containing deposits which may be present after measurement in dairy products such as milk, yogurt, cheese etc.; they contain pepsin in dilute hydrochloric acid. The membrane must never be treated with abrasive cleaners (Ata, Vim etc.).

- 1.5 If the electrode is used for pH measurement in non-aqueous media it should be soaked in water between tests in order to regenerate the soak layer of the glass membrane.

#### 2. Metal electrodes

- 2.1 The metal surface must be free from any contamination. Cleaning methods, apart from ultrasonics, grinding and polishing, may include the chemical agents indicated under II.1.4; in addition chrome-sulfuric acid to remove fatty layers.

#### 3. Reference electrodes

- 3.1 Only 3 M KCl solution should be used for storage or conditioning (or the appropriate electrolyte in the case of double-chamber electrodes) in order to avoid any reduction in salt concentration of the reference cell.
- 3.2 Reference electrodes with liquid electrolyte (identified by the blue piece of tubing closing the filler opening) must always contain sufficient liquid. During calibration and measurement the filler opening is kept open to ensure pressure equilibration; the electrolyte is topped up if necessary. In the case of measurements without external pressure the liquid column inside the electrode should be a few centimetres higher than the surface of the test solution in order to generate a positive hydrostatic pressure which ensures a continuous flow of KCl solution through the diaphragm outwards into the test solution. This ensures continuous self-cleaning of the diaphragm and in addition prevents diffusion of possible electrode poisons into the interior of the electrode which would render the measuring system unusable.

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3.3 In reference electrodes with gel reference electrolyte (identified by a natural-colour tubing piece closing the filler opening) no topping up is required. Such electrodes are maintenance-free. The tubing piece must not be removed.

3.4 With dry storage, potassium chloride crystallises out on the outside of the diaphragm. The electrodes must be rinsed thoroughly with water before use, or in case of a blocked diaphragm (indicated by drifting of the reading) stored for some time in 3 M KCl solution. In particularly difficult cases it is advisable to warm up the electrode in a water bath at a maximum temperature of 80°C. When an electrode is out of use for a longer period, the protective cap filled with KCl solution must be placed back over the electrode.

3.5 For cleaning dirty diaphragms the same methods as for glass electrodes can be used. In case of hard deposits there is in addition the possibility of cleaning the diaphragm surface with emery paper or by grinding. The pH glass membrane must not be scratched during cleaning.

#### 4. pH or metal combination electrodes

4.1 The details described for the appropriate single electrode (measurement or reference electrode) apply also to combination electrodes.

#### 5. Warranty

5.1 The usefulness of measurement and reference electrodes depends on their maintenance and on the conditions of use. The period until an electrode has to be replaced can vary between a few weeks and several years depending on the particular application.

5.2 In case of any complaint please return the faulty electrode which must still carry the serial number (label on the electrode head). In addition we require details of the fault.

5.3 1 year warranty. A quality product deserves confidence. For our electrodes we offer a 12 months' warranty on material and manufacturing faults.

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### III. Calibration and measurement

#### 1. Redox potential

The redox potential is indicated directly in millivolt; the meter must therefore be set to voltage. The electrode has to be connected up in accordance with the Operating Instructions of the instrument manufacturer. Special redox buffer solutions are available for testing the function of the measuring system.

#### 2. pH

A calibration must always be performed when using a new pH electrode system for the first time. In case of electrodes in continuous use it is advisable to check the function at regular intervals combined with recalibration.

#### 3. Calibration

The calibration instructions of the instrument manufacturer have to be observed. Where no instructions are available the following procedure has to be followed: Since electrodes are subject to certain individual tolerances and in addition exhibit drift it is advisable to perform a 2-point calibration. This requires two standard buffer solutions. During calibration it is important that before and after each buffer change the electrode is rinsed either in de-ionised water (followed by dabbing with a soft paper tissue) or in the subsequent buffer solution.

3.1 The electrode system is connected to the pH meter in accordance with the instrument Operating Instructions. The electrodes must be immersed in the solution so that both glass membrane and diaphragm are surrounded by solution.

3.2 The electrode system is immersed in a buffer solution whose pH is as close as possible to the electrode zero (pH<sub>0</sub> = 7.0). The temperature control on the meter is set to the solution temperature. For very accurate measurement (±0.01 pH) a heating bath must be used to bring the solution to the subsequent operating temperature of the electrode. Using the zero control the indication is set to the pH of the buffer solution.

3.3 The electrode system is now immersed in a second buffer solution at the same temperature; its pH (either

acidic, e.g. pH 4.00, or alkaline, e.g. pH 9.00) should be as near as possible to that of the solution in which the electrode will later be used. The pH of this second buffer solution is set on the scale using the slope potentiometer.

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#### 4. Measurement

To measure the pH of a solution the pH electrode system is immersed in the solution, the solution temperature is set on the meter, and the pH is read on the scale.

### IV. Faults

All electrodes are carefully tested before they leave the factory. Failure in practical use cannot however be avoided completely. These failures may be due either to the electrodes or to other components of the installation. Unsatisfactory treatment of electrodes may also lead to failure.

The checking of returned electrodes is a very costly procedure. It is therefore absolutely essential that you provide full details of the type of fault.

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