DinoEye eyepiece cameras HR 5 megapixel









 Cell structure viewed with eyepiece camera
Tilia stem viewed with eyepiece camera

HIGH RESOLUTION EYEPIECE CAMERA'S





Can be mounted on an endoscope when combined with a C-Mount coupler

ACCESSORY



The C-mount coupler is used to connect a DinoEye eyepiece camera to an endoscopic camera. The C-mount adapter has a standard 25mm diameter male thread.

For proper measurement results, it is recommended to always perform a calibration.